



MODEL 200

Dimpling Grinder

An easy-to-use, state-of-the-art, mechanical thinning instrument designed for the reproducible preparation of high-quality electron microscopy specimens.

Model 200 Dimpling Grinder specifications

Grinding control	Automated grinding rate control Grinding force adjustable with micrometer counterweight system Independent control of grinding wheel and specimen rotation speeds
Specimen stage	Precise specimen stage rotation Magnetic mount allows easy specimen positioning
User interface	All program inputs via front panel keypad Specimen thickness indicated on a liquid crystal display
Specimen illumination	Specimen observation in either transmitted or reflected light
Enclosure	Weight: 8.2 kg [18 lb.] Size: 208 mm [8.2 in.] width x 165 mm [6.5 in.] height x 343 mm [13.5 in.] depth
Power requirements	110/220 VAC, 50/60 Hz, 375 W
Warranty	One year



E.A. Fischione Instruments, Inc.
9003 Corporate Circle
Export, PA 15632 USA
Tel: +1 724.325.5444
Fax: +1 724.325.5443
info@fischione.com
www.fischione.com

©2024 E.A. Fischione Instruments, Inc. All rights reserved.
Document Number SP200 Revision 02 09/2024